


# ECLIPSE

AMNION MEMBRANE TISSUE

## SAFE. NATURAL. HEALING

### ECLIPSE ADVANTAGES

- Eclipse is air-dried not freeze-dried
- Eclipse can be kept at room temperatures, eliminating the need for freezer
- Suitable for immediate use off the shelf
- Hydrates rapidly
- Adheres to the ocular surface
- Sterilized with gamma radiation to SAL  $10^{-6}$  in accordance with ISO 11137
- Textured surface facilitates ease of handling and placement
- Symmetric graft allows either side to be placed in contact with ocular surface
- Can be placed in either direction (multi-directional graft)
- Has a micron thickness of 30-35 Microns
- Tackier stromal surface for better adherence
- Only 9 units of KiloGray used to sterilize the grafts
- Low doses of radiation retain more growth factors
- New foil packaging offers a 5-year shelf life
- Available in 5 mm, 8mm, 10 mm, 12 mm, 14 mm, 16 mm diameters
- More growth factors and cytokines to survive the sterilization process



**OPHTHALOGIX**  
an Evologics Company

[www.ophthalogix.com](http://www.ophthalogix.com)

OPT-AUR02-REV01

# AURORA AMT

INNOVATIVE, INTELLIGENT AMT POWERED BY OPHTHALOGIX



## AURORA SIZES



12 mm with 5mm aperture  
EAS-90012 - Aurora Single Layer 12mm  
EAD-90012 - Aurora Dual Layer 12mm



10 mm with 5mm aperture  
EAS-90010 - Aurora Single Layer 10mm  
EAD-90010 - Aurora Dual Layer 10mm

## AURORA, INNOVATIVE AMT

Better for Your Patients, Better for Your Practice – Healing Eyes, Naturally

At Ophthalmogix, we are committed to pushing the boundaries of medical innovation to provide your patients with the best possible outcomes.

Embrace the future of eye care with Aurora – where medical innovation leads to better patient outcomes – naturally. Learn more about how Aurora can improve patient quality-of-life.

Aurora's groundbreaking feature is a strategically placed hole, or aperture, in the center of the tissue. The aperture spares patients the reduced vision post-procedure associated with 1st generation grafts. This aperture permits maximal oxygen flow to the central cornea, supporting corneal wellbeing and eliminating the risk of central hypoxic keratitis and its sequelae.

## FORMULA/5™ Processing Technique

The FORMULA/5™ processing method is designed to retain natural growth factors and other regulatory proteins to support the full healing cascade.  
1,2,3,4,5

- Ensures minimal manipulation
- Uses no harsh chemicals
- Provides the most gentle drying process in the industry
- Engineered to meet stringent FDA regulations and AATB standards
- Room-Temperature Storage
- Quick Rehydration
- Easy Application

1. Partial list of regulatory proteins. 2. Data on File. Ophthalmogix can provide information upon request from C5 Biomedical, manufacturer of Aurora. 3. Koizumi N, et al. Current Eye Research. 2000; 20(3): 173-77. 4. Koob TJ, et al. J Biomed Mater Res Part B 2015; 103B: 1133-40. 5. Li JF, et al. HGF accelerates wound healing by promoting the dedifferentiation of epidermal cells through 1-integrin/ILK pathway. Biomed Res Int. 2013;2013:470418. doi: 10.1155/2013/470418. Epub 2014 Jan 15. PMID: 24490163; PMCID: PMC3899705.



**OPHTHALOGIX**  
an Evologics Company

[www.ophthalmogix.com](http://www.ophthalmogix.com)  
888-77-AMNION

